

Liko Universal Sling (Mod. 000)



- Liko UniversalSling is a basic model which is designed to adapt to the patient without individual adjustments and is one of our most used slings
- The UniversalSling provides an upright sitting posture and supports the entire back up to the neck
- The patient's arms are held inside the sling
- A properly fitted and carefully applied UniversalSling will give the patient an excellent sense of safety, as well as a high degree of comfort

Warranty Information

- 1 Year



✓ Contract approved

NDIS

Item Code:

05_123603164_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS394775	Universal Sling - Reinforced Leg Support - Polyester - Medium	200 kg	Polyester
LSS394774	Universal Sling - Reinforced Leg Support - Polyester - LARGE	200 kg	Polyester
LSS394773	Universal Sling - Reinforced Leg Support - Polyester - XLARGE	300 kg	Polyester