

Paragolfer



ottobock.

- The Paragolfer is a powerful Sit-to-Stand powerchair for use in Sports, primarily but not limited to golf.
- The wide base of 90cm ensures you have a stable base for that perfect golf swing, while the broad drive and trailing wheels help to protect the green.

✓ Contract approved

NDIS

Item Code:

05_122203291_0112_1_2

Wheelchair Powered –

Composite:05_419_0105_1_2

Click [here](#) for more information

Code	Display Name	Total Length	Seat Width	Total Width	Safe Working Load	Maximum Speed
PWC654230	Paragolfer	149 cm	46 cm	90 cm	140 kg	10 km/h