

CP200 Portable Ceiling Hoist



PRISM MEDICAL

- Hoist can be reversed to act as a portable or fixed unit providing flexibility of use
- Slimline, low profile, padded carry bar that has six loop attachments as standard (four-point hook optional). Different width carry bars also available
- Low profile hoist unit is unobtrusive and provides greater lifting height
- All lifting and traversing operations have soft start/ stop to ensure smooth transfers for the patient
- A host of safety features
- Smooth running tape with large lifting range
- State of the art display with lift counter
- Ergonomically designed handset with hold to run clearly marked buttons and directional arrows requires minimum pressure to operate. Pneumatically controlled, making it safe for use in wet conditions
- Fire retardant cover completely encapsulates hoist
- Neat in track transformer, IP 67 rated
- Coloured direction indicators on hoist unit match hand control to ensure correct operation first time
- Emergency stop and manual emergency window

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Warranty Information

- 5 Years

Code	Display Name	Safe Working Load
------	--------------	-------------------

LIC350105 Prism CP200 Ceiling Hoist 200 kg

Code	Display Name	Safe Working Load
------	--------------	-------------------