

PAWS - Accessories



- **STORAGE BAG** Secures firmly with velcro fasteners. Zip top and side pockets keep valuables safe and essentials close to hand.
- **STORAGE BASKET** Easy clip on and off basket that allows you to load up for your trip. With sturdy heavy gauge steel rod construction, it is dip coated to prevent corrosion.
- **PHONE HOLDER** Allows phone to be safely stored and visible when in use. Facilitating hands free functions such as navigation.
- **SMART TYRE INFLATOR & PHONE CHARGER** The mini compressor with rechargeable battery can automatically inflate tyres to specific pressures and offers USB charging for other devices.
- **QUICK SWAP BATTERIES**
 - **Lithium Ion Battery Pack – 11.6AH – WCA872220** Secured under lock and key, the powerful lithium ion battery packs can be conveniently swapped avoiding the need to wait for charging. Full charge time is 5 hours and each has an inbuilt charge indicator.
 - **Lithium Ion Flight Battery Pack – 5.6AH – WCA872230** Air travel is possible with the optional Flight battery. IATA approved for travel on passenger planes, this smaller battery pack meets the maximum power levels allowed by commercial carriers.



Contract approved

Support @ Home

Function: Wheelchair accessories

Category: Prescribed

Click [here](#) for more information

Code	Display Name
WCA872210	PAWS - Smart Tyre Inflator & Phone Charger
WCA872220	PAWS - Lithium Ion Battery Pack – 11.6AH

Code	Display Name
WCA872230	PAWS - Lithium Ion Flight Battery Pack – 5.6AH
WCA872260	PAWS - Phone Holder
WCA872270	PAWS - Storage Basket
WCA872280	PAWS - Storage Bag
WCA872200	Quad/Tetra Hand Control - Rehasense Paws