

TDX SP2



The Invacare TDX SP2 is a high performance, centre wheel drive, power wheelchair equipped with the advanced Ultra Low Maxx seating system and revolutionary LiNX smart technology. This power wheelchair is the perfect mix of functionality, customisation and style that adapts to the user's needs and desired lifestyle.

- Enhanced stability: The Stability Lock function ensures all six wheels stay firmly on the ground on uneven terrain
- LiNX Smart Technology: The smart control system enables programming and maintenance updates to be made wirelessly and in real time, plus it has straightforward access to statistical and diagnostic information, giving complete visibility to professionals
- Range of specialist controls
- Intelligent seating system: The Ultra Low Maxx seating system can combine up to 50° tilt and 170° of recline
- Intuitive driving experience
- Centre Wheel Drive (CWD) provides the smallest turning diameter and smallest overall footprint vs FWD and RWD for optimum maneuverability
- Find premium solutions for the most clinically complex needs with Motion Concepts Ultra Low Maxx (ULM) Power Positioning System: tilt, elevate and power legrests
- G-Trac Technology helps provide precise tracking for imprecise environments such as side slopes, uneven or soft ground, or even bunched up carpets or door jambs.
- Invacare LiNX Technology is an insight inspired control system with advanced technology that provides a superb driving experience for users and allows professionals to configure and tailor the TDX SP2 Power Wheelchair quickly and more intuitively
- SureStep Suspension and Stability Lock work in synergy to help maintain

✓ Contract approved

Support @ Home

Function: Powered wheelchairs

Category: Prescribed

Click [here](#) for more information

balance and provide stability when transitioning over obstacles

Code	Display Name
PWC652317	TDX SP2 (as per script)