

Pride Pathrider 130XL Pursuit Scooter



PRIDE®
MOBILITY
LIVE YOUR BEST®

- Large 33 cm (13") pneumatic tyres for excellent outdoor performance
- Full lighting package including directional signals and angle-adjustable, dual bulb lower headlight
- Seat post suspension for a smoother ride
- Easy drive tiller with wraparound handles
- Front basket standard
- Easily accessible tie down points (for transport of unoccupied scooter)
- Rearview mirror standard

Warranty Information

- 2 Years
- 6 months - Batteries

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Diverse motorised vehicles

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_122303123_0105_1_2

Scooter –

Composite:05_418_0105_1_2

Click [here](#) for more information

Compliance

AS3695.2:2013

Code	Display Name	Total Length	Total Width	Safe Working Load	Maximum Speed	Range	Colour
SCT670571	Pride Pathrider 130XL PURSUIT - 4 Wheel - Red	1360 mm	635 mm	159 kg	12.4 km/h	25km	Red

Code	Display Name	Total Length	Total Width	Safe Working Load	Maximum Speed	Range	Colour
SCT670575	Pride Pathrider 130XL PURSUIT - 4 Wheel - Yellow	1360 mm	635 mm	159 kg	10 km/h	25 km	Yellow