

Liko Viking XL Lifter



- Ideal for bariatric lifts and transfers
- Universal TwinBar 670 Slingbar offers additional sling space for patients of size
- Informative control box displays number of lift cycles, battery level, lift overload and service indicator
- Optional quick-release hook system which means caregivers can easily and rapidly change a slingbar or install a scale
- Optional armrests help support patients during ambulation as well as making it easier for caregivers to manoeuvre the lift
- Three variable height settings accommodate a variety of different lifting situations
- Intuitive and ergonomically-designed hand control incorporates indicators for low battery level and when the Safe Working Load of the lift is exceeded
- An extensive assortment of accessories is available including slings, slingbars and scales
- Lifting Height 540 - 1910 mm

Warranty Information

- Frames: 10 Years
- Electrical Parts: 2 Years

 Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_123603161_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material
LSS393170	Liko Viking Lifter - XLARGE	1420 mm	815 - 1130 mm	1500 - 2110 mm	300 kg	Aluminium

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material
LSS393220	Liko Viking XL - Leg Protector Set -	-	-	-	-	-