

JAY - Basic



JAY®

- Designed for the client who has symmetrical postures needing minimal postural support
- Durable pre-contoured foam provides lateral stability; encourages orthopaedic alignment; and helps redistribute pressure to increase sitting tolerance
- Most suitable for clients with low risk of skin breakdown
- All models include: Incontinence resistance cover

Warranty Information

- 1 Year Cushion
- 6 Months Cover

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_181006190_0105_1_2

Pressure Care Cushion –

Composite:05_429_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Total Height	Safe Working Load	Pressure Care Rating	Material
PTC567600	JAY Basic - 14 x 16in	2.5 in	150 kg	Low Risk	Foam
PTC568400	JAY Basic - 15 x 16in	2.5 in	150 kg	Low Risk	Foam

Code	Display Name	Total Height	Safe Working Load	Pressure Care Rating	Material
PTC567400	JAY Basic - 16 x 16in	2.5 in	150 kg	Low Risk	Foam
PTC567500	JAY Basic - 16 x 18in	2.5 in	150 kg	Low Risk	Foam
PTC568450	JAY Basic - 17 x 16in	2.5 in	150 kg	Low Risk	Foam
PTC566800	JAY Basic - 18 x 16in	2.5 in	150 kg	Low Risk	Foam
PTC566900	JAY Basic - 18 x 18in	2.5 in	150 kg	Low Risk	Foam
PTC567650	JAY Basic - 20 x 16in	2.5 in	150 kg	Low Risk	Foam
PTC567700	JAY Basic - 20 x 18in	2.5 in	150 kg	Low Risk	Foam
PTC567750	JAY Basic - 22 x 18in	2.5 in	150 kg	Low Risk	Foam
PTC566950	JAY Basic - 22 x 20in	2.5 in	150 kg	Low Risk	Foam