

Malvestio Idea 2 Treatment Chair



- Premium therapy chair that effortlessly meets the needs of oncology departments and the demands of infusion and transfusion medicine.
- Gas-assisted, height and laterally adjustable swing away armrests to suit individual user needs.
- Electric adjustment allows for infinite positioning of the backrest, legrest, seat height, seat tilt and Trendelenburg (CPR shock position).
- Ideal for procedures such as chemotherapy, dialysis, blood donation, and blood transfusion



Code	Display Name	Total Length	Total Width	Safe Working Load
CHT100710	Idea 2 Treatment Chair	1600 mm	840 mm	200 kg