

Liko Sabina II EE Lifter



- Conventional patient lifting to a limited extent, lifting in the sitting position
- Powerful lifting motor ensures the smooth raising and lowering of patients
- Base-width adjustment offers greater flexibility of use when lifting from a bed or chair, enhancing patient mobility
- The lower leg support and the foot rest can easily be removed for cleaning
- Adjustable height settings accommodate a variety of different patient sizes
- Sling bars are sold separately to give customer the choice of best one for their needs. We recommend either the Sabina II 350 SlingBar (LSS393270) or Sabina II Comfort Sling Bar (LSS393280)

Warranty Information

- Frames: 10 Years
- Electrical Parts: 2 Years



Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_123604167_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material
LSS393250	Liko Sabina II EE Lifter	1060 mm	695 - 1115 mm	900 - 1000 mm	200 kg	Steel
LSS393320	Liko Sabina Seat Strap SlingBar	-	-	-	-	-
LSS393330	Liko Sabina II Seat Strap SlingBar	-	-	-	-	-

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material
LSS393340	Liko Sabina Seat Strap	-	-	-	-	-
LSS393350	Liko Sabina Calf Strap	-	-	-	-	-
LSS393270	Liko Sabina II 350 SlingBar	-	-	-	-	-
LSS393280	Liko Sabina II Comfort SlingBar	-	-	-	-	-
LSS393290	Liko Sabina Side Support	-	-	-	-	-
LSS393300	Liko Sabina II Side Support	-	-	-	-	-
LSS393310	Liko Sabina II Heel Support	-	-	-	-	-