

Liko Hygiene Sling (Mod. 40 & 41)



- The Liko Hygiene Sling is designed in order to assist with dressing and undressing when lifting to and from the toilet, these are lifting situations that are often difficult for caregivers
- Thanks to the large opening around the buttocks, dressing and undressing can be performed during the lifting procedure itself
- The Hygiene Sling is easy to apply even if the patient is placed tightly in the wheelchair. The patients arms should always be on the outside of the Hygiene Sling to ensure safety
- Model 41 is equipped with a belt for added safety and security
- The Hygiene Sling usually works well for patients with long trousers and for patients wearing skirts and leggings

Warranty Information

- 1 Year



Contract approved

NDIS

Item Code:

05_123603163_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS394480	Hygine Sling Mod 41 - Belt and Reinforced Leg Support - Polyester - SMALL	200 kg	Polyester (Loops)
LSS394490	Hygine Sling Mod 41 - Belt and Reinforced Leg Support - Polyester - MEDIUM	200 kg	Polyester (Loops)
LSS394500	Hygine Sling Mod 41 - Belt and Reinforced Leg Support - Polyester - LARGE	200 kg	Polyester (Loops)
LSS394510	Hygine Sling Mod 41 - Belt and Reinforced Leg Support - Polyester - XLARGE	300 kg	Polyester (Loops)