

Q6 Edge 2.0 / iLevel



- 9.6 km/h high-speed motors standard
- Responsive handling Precise control
- Advanced performance
- Bold aesthetics
- Optional iLevel features Extra Stability Technology to ensure maximum safety while performing transfers, reach, and numerous other activities of daily living iLevel enables up to 10" of seat elevation while at walking speed (up to 5.6 km/h)
- iLevel is retrofittable on the Q6 Edge®
- 2.0 12 color options Available with Accu-Trac Advanced Tracking Technology
- Optional LED fender lights Optional 8" casters'



Contract approved

Support @ Home

Function: Powered wheelchairs

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_122306138_0105_1_2

Wheelchair Powered –

Composite:05_419_0105_1_2

Click [here](#) for more information

Code	Display Name	Total Length	Seat Width	Total Width	Safe Working Load	Maximum Speed
PWC652226	Q6 Edge 2.0 / iLevel	90.2 cm	30 - 56 cm	62 cm	136 kg	9.66 km/h