

Aspire Raphael CUSTOM Quattro Powerlift Recline Chair - Soft Fabric



- The soft touch weave and diverse colour palette makes the Mazza and Manisa ranges a highly practical and functional solution
- These fabrics are comforting, easy to clean and durable, yet still soft to touch
- Commercial quality fabric for long lasting durability
- Other custom order fabrics include [Lustrell Ranger](#) (Healthcare Vinyl), [Eastwood](#) (Micro Suede) and [La Casa](#) (Leather)
- Custom order lead time varies, approximately 5 - 6 weeks depending on fabric availability
- Custom orders are subject to terms and conditions, please see page two on the order form below
- All colour swatches shown may not be a true representation of the actual product and should be used as a guide only
- Ready available Raphael Powerlift Recline Chair [options here](#)
- All colour swatches shown may not be a true representation of the actual product and should be used as a guide only
- **Raphael Powerlift Recline Chair Features & Benefits:**
 - Infinite positioning can be achieved with the Raphael Powerlift Recline Chair thanks to its 4-motor engineering
 - A true centre-piece in the modern home
 - Backlit hand control with an integrated USB charging point
 - Independently adjustable headrest, back, leg & lumbar support functions
 - Integrated lift & recline functions for easy sit-to-stand & maximum comfort
 - Intuitive headrest support, ideal when watching TV when reclined
 - Available in four sizes: Small, Medium, Large and MAXI

✓ Contracted DVA product



Warranty Information

- 2 Years

Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
CHP227290FA	Raphael CUSTOM Quattro Powerlift Recline Chair - Soft Touch Fabric - SMALL	740 mm	1060 mm	130 kg	480 mm	480 mm	470 mm	470 mm	630 mm
CHP227300FA	Raphael CUSTOM Quattro Powerlift Recline Chair - Soft Touch Fabric - MEDIUM	780 mm	1090 mm	130 kg	500 mm	520 mm	470 mm	470 mm	680 mm
CHP227310FA	Raphael CUSTOM Quattro Powerlift Recline Chair - Soft Touch Fabric - LARGE	810 mm	1080 mm	130 kg	530 mm	550 mm	500 mm	500 mm	700 mm
CHP227320FA	Raphael CUSTOM Quattro Powerlift Recline Chair - Soft Touch Fabric - MAXI	930 mm	1050 mm	180 kg	540 mm	650 mm	510 mm	510 mm	650 mm