

Aspire 460mm Weight Bearing Platform Shower Commode



- 8° backrest angle ensures patient comfort
- Improved over toilet clearance
- Integrated access trough improves hygiene whilst maintaining toilet bowl clearance height of 465 mm
- 4 x Tente total locking castors provide brake access from all angles
- Safety arms prevent falls and swing back out of the way for transfers
- Open armrests style accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility
- Easy-slide weight bearing platform footrest slides completely away for transfers
- Overall Height: 995 mm
- Overall Width: 535 mm
- Seat to Backrest: 305 mm
- Armrest Height: 225 mm
- Anterior Seat Height: 590 mm
- Posterior Seat Height: 560 mm
- ***Seat sold separately**
- Note: Seat is supplied automatically with any rental equipment

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_091203053_0103_1_2

Mobile Shower Commode –

Composite:05_413_0103_1_2

Click [here](#) for more information

Prescribing Therapists: Click here to download our
[Clinical Assistance Guide](#)

Warranty Information

- Corrosion: 5 Years
- Components: 1 Year

Compliance

AS/NZS 3973:2009

Code	Display Name	Seat Width	Safe Working Load
BTC066000	Attendant Propelled Commode - Weight Bearing Platform - 460 mm Wide	460 mm	200 kg