

Aspire Socialite Folding Wheelchair - Attendant Propelled



- The Aspire Socialite is an ultra light weight transport chair perfect for short trips and outings with friends.
- The generous 115kg safe working load caters to a broad range of users.
- Padded and breathable seat and back provide optimum comfort and postural support.
- Cut away arms improve table access and the offset uprights create hip space for coats and warmer clothing.
- Sleek European design.
- Stylish contemporary colours.
- One of the lightest transport chairs available.
- Perfect for outings and appointments.
- Attendant operated hand brakes with parking lock.
- Simple and intuitive folding for easy transportation.
- Removable leg rests and folding back for compact storage.
- Contoured seating improves comfort and postural support.
- Rear storage pocket for essentials.
- Seat height anterior: 475 - 500 mm.
- Seat height posterior: 440 mm.
- 5 Years warranty - Frame
- 1 Year - Fabric & Components

✓ Contract approved

Support @ Home

Function: Manual wheelchairs

Category: Prescribed

Click [here](#) for more information

Compliance

Australian Standard AS/NZS

3695.1:2011

Code	Display Name	Total Width	Safe Working Load	Unit Weight	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Frame Weight	Backrest Height
MWS449800	Aspire SOCIALITE Folding Wheelchair -	595 mm	115 kg	8.9 kg	Silver	430 mm	450 mm	440 mm	440 mm	8.1 kg	410 mm

Code	Display Name	Total Width	Safe Working Load	Unit Weight	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Frame Weight	Backrest Height
	Attendant Propelled - Silver										
MWS449810	Aspire SOCIALITE Folding Wheelchair - Attendant Propelled - Red	595 mm	115 kg	8.9 kg	Red	430 mm	450 mm	440 mm	440 mm	8.1 kg	410 mm