

Aspire Rembrandt Quattro Lift Recline Chair



The Rembrandt Quattro offers a unique seating experience. With independent control of back, leg, lumbar and headrest functionality, the user has infinite positioning options with the press of a button. Genuine shiatsu massage coupled with optional heating gives users a truly unique and luxurious comfort experience. Zero Gravity tilt-in-space positioning creates a feeling of weightlessness and promotes pressure redistribution.

- **Shiatsu Roller Massage** in backrest provides therapeutic and stress relief for the user
- **Memory Function** allows for a preset favourite position to be set
- **Independent Adjustable Headrest** improves comfort whilst watching TV or reading a book
- **Integrated Lift and Recline** functions for easy sit-to-stand
- **Extending Legrest** supports more of the lower legs for enhanced comfort
- **Heating Function** within the seat, lumbar and backrest provides warmth and comfort during the colder months
- **Cup Holder** built into the armrest for easy drink storage
- **Zero Gravity** tilt-in-space positioning creates a feeling of weightlessness and promotes pressure redistribution
- **Wireless Device Charging** built into the armrest for convenience
- **Lithium Ion Battery Backup** ensures chair can be used in the event of a power outage or if moved between rooms during the day
- **Backlit Hand Control** with integrated USB charging point

Warranty Information

- 5 Years - Frame
- 2 Years - Components



Contract approved

Support @ Home

Function: Sitting furniture

Category: Under advice

Click [here](#) for more information

Code	Display Name	Total Width	Total Height	Safe Working Load	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height	User Height
CHP227550	Aspire Rembrandt Quattro Lift Recline Chair	1015 mm	1095 mm	159 Kg	Canyon Walnut	500 mm	560 mm	495 mm	495 mm	600 mm	1065 mm