

## Mygo Stander



**LECKEY®**

- The Mygo Stander provides 3-in-1 prone, supine and upright standing positions
- Comes with height adjustable chest pad that is ergonomically designed to imitate the shape of the rib cage
- Includes hip support, chest support, 4 castors with brakes and height adjustable footplates and knee supports
- Chest Width - PED479335 = 160 - 280 mm, PED479336 = 200 - 320 mm
- Hip Width - PED479335 = 220 - 330 mm, PED479336 = 280 - 390 mm
- Height of User Range - PED479335 = 1080 - 1450 mm, PED479336 = 1350 - 1750 mm

### Warranty Information

- 2 Years



**NDIS**

**Item Code:**

05\_053603010\_0103\_1\_2

Click [style='color: #6c0577;'here](#)

for more information

| Code      | Display Name                            | Total Width | Safe Working Load |
|-----------|-----------------------------------------|-------------|-------------------|
| PED479335 | Mygo Stander, Size 1, Ages 4 - 10 years | 630 mm      | 50 kg             |
| PED479336 | Mygo Stander, Size 2, Ages 8 - 14 years | 650 mm      | 60 kg             |