

Aspire Socialite Folding Wheelchair - Self Propelled



- The Aspire Socialite is an ultra light weight transport chair perfect for short trips and outings with friends.
- The generous 115kg safe working load caters to a broad range of users.
- Padded and breathable seat and back provide optimum comfort and postural support.
- Cut away arms improve table access and the offset uprights create hip space for coats and warmer clothing.
- 22in rear wheels for user self-propulsion.
- Sleek European design.
- Stylish contemporary colours.
- One of the lightest transport chairs available.
- Perfect for outings and appointments.
- Attendant operated hand brakes with parking lock.
- Simple and intuitive folding for easy transportation.
- Removable leg rests and folding back for compact storage.
- Contoured seating improves comfort and postural support.
- Rear storage pocket for essentials.
- Folded Dimensions (DxWxH) 650x380x740 mm
- Seat height anterior: 475 - 500 mm.
- Seat height posterior: 440 mm.

Warranty Information

- 5 Years - Frame
- 1 Year - Fabric & Components



Contract approved

Support @ Home

Function: Manual wheelchairs

Category: Prescribed

Click [here](#) for more information

Standards

Australian Standard AS/NZS

3695.1:2011

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Depth	Seat Width	Backrest Height
MWS449820	Aspire SOCIALITE Folding Wheelchair - Self Propelled - Silver	1030 mm	675 mm	900 mm	115 kg	13 kg	Silver	430 mm	450 mm	410 mm
MWS449830	Aspire SOCIALITE Folding Wheelchair - Self Propelled - Red	1030 mm	675 mm	900 mm	115 kg	13 kg	Red	430 mm	450 mm	410 mm