

Aspire Autoclavable Bed Pan



Our Autoclavable Bed Pans are made from premium polypropylene. They will last for many years in any healthcare environment subject to the rigours of cleaning, sanitation and infection control requirements.

- Heavy duty and reusable
- Designed for use in healthcare environments such as surgeries, acute care, general ward, medical centres, and aged care facilities
- Autoclavable polypropylene that can be cleaned using autoclaves and sanitisers
- Pastel green colour for easy identification

Warranty Information

- 1 Year



✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Collecting urine and faeces

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_098800044_0103_1_2

Toilet Aids – Composite:

05_414_0103_1_2

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Code	Display Name	Total Length	Capacity
BTT151800	Bed Pan	360 mm	1000 ml

