

Pride Rangerider Scooter



PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

- Rugged and built to last, the Rangerider is everything you want in an outdoor scooter!
- Equipped with front and rear suspension, this robust, rear wheel drive scooter is designed for cruising on walking trails or driving around your property
- Enjoy total visibility with complete LED lighting
- With a user-friendly console and dual hydraulic brakes for added safety
- Limited recline, black vinyl, high-back seat with lumbar and pillow-top headrest, flip-up height, and angle-adjustable armrests
- Scooter Frame: 2 Years
- Motors/Axles/Brakes: 2 years
- Battery: 6 Months

✓ Contract approved

Compliance

AS3695.2:2013

Code	Display Name	Total Length	Total Width	Safe Working Load	Colour	Maximum Speed	Range
SCT670581	Pride Rangerider Scooter - 4 Wheel	1560 mm	830 mm	159 kg	Grey	10 km/h	30 - 36 km