

Revolution Chair



- The Revolution Chair incorporates a unique, patented “swivel-slide” mechanism
- Swivel – to allow easy access to the chair seat
- Slide – to bring you close to the table edge
- At all times, the chair legs remain fixed and stable
- Sturdy steel frame, waterfall polyurethane arms
- Bleach-cleanable commercial grade vinyl seat and back

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Sitting furniture

Category: Under advice

Click [here](#) for more information

Code	Display Name	Total Width	Safe Working Load	Colour	Seat Depth	Seat Width	Seat To Floor Height
CHP196832	Revolution Chair - Height Adjustable - Onyx	505 mm	120 kg	Onyx	485 mm	450 mm	450 - 530 mm
CHP196834	Revolution Chair - Height Adjustable - Sand	505 mm	120 kg	Sandal Fabric	485 mm	450 mm	450 - 530 mm
CHP196830	Revolution Chair - Height Adjustable (Custom - Specify Colour)	505 mm	120 kg	-	485 mm	450 mm	450 - 530 mm
CHP196820	Revolution Chair - (Custom - Specify Colour)	505 mm	120 kg	-	485 mm	450 mm	480 mm
CHP196840	Revolution Chair - Bariatric (Custom - Specify Colour)	800 mm	250 kg	-	480 mm	660 mm	500 mm