

Aspire Vogue Adventure Walker



- **Heavy-duty** frame and wider seat provide enhanced support for users up to 200kg.
- **European Design** offers personalised height adjustment, promoting correct posture.
- **Large Multi-Terrain Wheels** for durability and stability on and off road. (Front Wheels 11", Back Wheels 9")
- **Ergonomic Handbrakes** lock for improved user safety and stability.
- **Compact Turning Circle** facilitates improved mobility in tight spaces.
- **Simple Folding Mechanism** enables easy transportation and compact storage (245 mm folded width, perfect for the car).
- **Detachable Storage Bag** carries up to 5kg of essentials.
- **Bonus** walking cane holder provides increased flexibility for your mobility needs.
- *Unit weight includes accessories, bag and backstrap
- **Frame weight is minimal functional weight

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Walking, manipulated by both arms

Category: Prescribed

Click [here](#) for more information

Compliance

ISO 11199-2:2005

Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Depth	Seat Width	Seat To Floor Height	Frame Weight
WAF705450	Aspire Vogue Adventure Walker	770 mm	695 mm	900 - 1035 mm	200 kg	10 kg	Black	230 mm	540 mm	530 mm	8.8 kg

