

Aspire ComfiMotion Memory Foam Mattress



- The ComfiMotion Memory Foam mattress contours to the user's body, providing personalised comfort and transitional pressure relief
- Integrated side-walls increase roll resistance at the bed edge and assist users, partners or carers with extra rigidity to ease bed entry and exit
- Contoured Immersion: Specially layered memory foam provides contoured immersion for support and pressure relief
- Superior breathability: Up to 8 x more breathable than conventional memory foams
- Anti-microbial: Environmentally safe treatment is added to the foam which reduces odours and kills 99.9% of bacteria
- Eco-friendly GECA Foam is produced using cleaner chemicals and production technology

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Compliance

BS7177:2008

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
BEM469910	Aspire ComfiMotion Memory Foam Mattress - Portobello Grey (sides) - Long Single	2030 mm	920 mm	200 mm	180 kg
BEM469930	Aspire ComfiMotion Memory Foam Mattress Pressure Care - Long Single	2030 mm	920 mm	200 mm	180 kg
BEM469920	Aspire ComfiMotion Memory Foam Mattress- Portobello Grey (sides) - King Single	2030 mm	1070 mm	200 mm	180 kg
BEM469940	Aspire ComfiMotion Memory Foam Mattress Pressure Care - King Single	2030 mm	1070 mm	200 mm	180 kg
BEM469922	Aspire ComfiMotion Memory Foam Mattress - Portobello Grey (sides) - Long Double	2030 mm	1370 mm	200 mm	360 kg
BEM469942	Aspire ComfiMotion Memory Foam Mattress Pressure Care - Long Double	2030 mm	1370 mm	200 mm	360 kg
BEM469924	Aspire ComfiMotion Memory Foam Mattress - Portobello Grey (sides) - Queen	2030 mm	1520 mm	250 mm	360 kg

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
BEM469944	Aspire ComfiMotion Memory Foam Mattress Pressure - Queen	2030 mm	1520 mm	200 mm	360 kg