


Aspire Slide Sheet



- Nylon low friction material assists in sliding and transferring residents
- Folded-over and placed under resident. The top sheet then slides over the bottom sheet during the transfer process
- Warranty: 12 months

Warranty Information

- 1 Year

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Changing body position

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_123103151_0105_1_2

Minor Transfer –

Composite:05_430_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width
TSE673070	Slide Sheet - Small	1000 mm	1450 mm
TSE673080	Slide Sheet - Large	2000 mm	1450 mm

