

B- Series



ottobock.

- The B-series power wheelchair has a modular design that offers almost boundless possibilities for assembling custom-made devices for complex requirements.
- With front-wheel and rear-wheel drive, the specific driving characteristics and motor powers of two kinds of drive in two frame widths are at your disposal.
- A high degree of flexibility is assured by the new VAS (variably adjustable seat) seat frame solution with Ergo-Joint back angle adjustment as well as innovative special and environmental controls, among other features.
- Optimised for more flexibility
- Back angle 0 to 30 degrees
- Seat Tilt 45 degrees power
- Seat inclination -3 / 0 / 3 / 6 (degrees)
- Max charge time of 12 hours

Warranty Information

- 2 Years



Contract approved

Support @ Home

Function: Powered wheelchairs

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_122306125_0105_1_2

Wheelchair Powered –

Composite:05_419_0105_1_2

Click [here](#) for more information

Code	Display Name	Seat Width	Total Width	Safe Working Load	Maximum Speed	Range
PWC652300	B-Series	340 - 560 mm	340 - 560 mm	160 kg	6 / 7.2 / 10 km/h	35 km