

## Pride Pathrider 140 XL Scooter



- Top of the range sleek and sporty powered mobility scooter
- Large 13" pneumatic tyres to handle rugged outdoor terrain
- Luxurious suspension seat with arm rests, head rest and lumbar support for stability
- Front and rear suspension and sealed hydraulic brakes
- Automatic full lighting system including directional adjustment. Long lasting high intensity LED kerb light, headlight, tail light and indicators
- Back lit battery gauge
- Wrap around easy drive tiller
- Tiller mounted battery charger port for off board charging
- Unit Weight: 118 kg

### Warranty Information

- 2 Years
- 1 Year Batteries

✓ Contracted DVA product

✓ Contract approved

### NDIS

#### Item Code:

05\_122303122\_0105\_1\_2

#### Scooter –

**Composite:**05\_418\_0105\_1\_2

Click [here](#) for more information

### Compliance

AS3695.2:2013

| Code      | Display Name                                | Total Length | Total Width | Safe Working Load | Maximum Speed | Range | Colour |
|-----------|---------------------------------------------|--------------|-------------|-------------------|---------------|-------|--------|
| SCT670578 | Pathrider 140 XL Scooter - Yellow - 4 Wheel | 1428 mm      | 692 mm      | 181 kg            | 10 km/h       | 35 km | Yellow |
| SCT670576 | Pathrider 140 XL Scooter - Red - 4 Wheel    | 1428 mm      | 692 mm      | 181 kg            | 10 km/h       | 35 km | Red    |
| SCT670577 | Pathrider 140 XL Scooter - Silver - 4 Wheel | 1428 mm      | 692 mm      | 181 kg            | 10 km/h       | 35 km | Silver |