

Aspire Go Turner



The Aspire Go Turner is a proven sit-to-sit transfer aid for quick and easy manual handling. This compact folding patient mover has a very narrow footprint, making it the perfect solution for transferring patients in tight spaces. The pull apart and folding design reduces storage space requirements and helps to transport the product if needed.

- Ergonomic handlebar
- Sling retention clips - Large range of sling belts (S to XL) for improved support and positioning
- Angle and height adjustable knee pads offer comfort during transport
- Central locking footbrake
- Spacious footplate
- 6 x wheels for improved stability

✓ Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Waist/Hip Circumference
LSS390585	Aspire Go Turner - Sit-to-Sit Transfer Aid	685 mm	570 mm	1205 mm	150 kg	18.9 kg	-
LSS391132	Aspire Go Turner Standing Belt - Small	-	-	-	150 kg	-	780 - 980 mm
LSS391133	Aspire Go Turner Standing Belt - Medium	-	-	-	150 kg	-	930 - 1180 mm
LSS391134	Aspire Go Turner Standing Belt - Large	-	-	-	150 kg	-	1130 - 1380 mm
LSS391135	Aspire Go Turner Standing Belt - Xlarge	-	-	-	150 kg	-	1330 - 1530 mm