

## Zippie Iris



- The Zippie IRIS® rotation-in-space wheelchair offers a superior combination of positioning, maneuverability, and increased portability with the XLOCK® folding crossbrace option.
- Cutting-edge rotation-in-space technology designed just for kids
- This tilting wheelchair uses patented rotation-in-space technology to rotate the seat frame around the user's center of gravity.
- This rotational system provides the shortest possible wheelbase without compromising stability to ensure a smooth and controlled rotation.
- The Zippie IRIS also comes with the True Fit Program, a complete growth package for free.

 Contract approved

### Support @ Home

**Function:** Manual wheelchairs

**Category:** Prescribed

Click [here](#) for more information

### NDIS

**Item Code:**

05\_122203115\_0105\_1\_2

**Wheelchair Manual –**

**Composite:**05\_417\_0105\_1\_2

Click [here](#) for more information

### Compliance

AS 3696-1:2008/AS

3696-3:2008/AS 3696-5:1989/AS

3696-8:1998/AS 3695:1992

| Code      | Display Name | Seat Depth | Seat Width | Seat To Floor Height | Safe Working Load | Unit Weight |
|-----------|--------------|------------|------------|----------------------|-------------------|-------------|
| MWC427500 | Zippie Iris  | 36 - 56 cm | 25 - 48 cm | 32 - 50 cm           | 125 kg            | 17 kg       |