

Extreme X8 - 4x4



magic
Mobility

- The Extreme X8 has a top speed of 10 km/h (6.2 mph)
- Articulating front frame with all terrain low pressure knobby tyres, which provides an exceptionally smooth and jolt-free ride
- Comes in five standard colours and three specialised prints
- Four independent motors to make kerbs and small steps easy to climb
- Splash resistant (not waterproof)



Contract approved

Support @ Home

Function: Powered wheelchairs

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Total Length	Total Width	Safe Working Load	Seat To Floor Height	Seat To Footrest	Maximum Speed
PWC655000	Extreme X8 - 4x4 (Base Price only)	1030 mm	710 mm	182 kg	470 mm	470 mm	10 kph