

PAWS Power Assisted Wheelchair Systems - Cruiser



- A great solution for a variety of mobility purposes. Whether it is around town or along a forest trail, this model is for you.
- The big 16" drive wheel eats up the bumps for a really smooth ride.
- Get your PAWS on the track and easily overcome the obstacles in the terrain.
- Just do more, further, faster and feel the freedom with PAWS.
- Optional Tetra Controls allow you to drive easily, even with diminished upper extremity strength and control.
- With PAWS your manual wheelchair becomes an exciting vehicle, allowing you to be independent when getting to work, going shopping, or on a trip.
- It combines ergonomics and functionality with good looks.
- Every PAWS creates its own story of freedom – meet family and friends, participate in more activities together or go for a solo adventure.
- When choosing PAWS you can customise the set up with different configurations of clamps, steering, wheels, and accessories for optimal use and fit – especially beneficial for quadriplec/tetraplegic users.
- With Auto Lift, you don't need to use upper body strength to connect PAWS with the wheelchair. Auto Lift will automatically connect you to driving position.
- PAWS comes with a variety of speed modes (*maximum speed is 10 km/h).
- PAWS has an IP54 waterproof rating so is not afraid of rain, puddles and moisture. Go on, get PAWS out!
- The Lithium Ion battery only needs up to 5 hours to be fully charged.
- On a full charge you can reach a range of up to 40 km (depending on terrain, driver weight and model chosen).
- An IATA compliant "flying" battery is also available for when you take your PAWS on a flight.
- With Automatic Easy Clamps just roll your wheelchair into the clamps and press the button for automatic clamping and lifting to the driving position.

✓ Contract approved

From Start to go in 12 seconds!

Code	Display Name
WCA872120	PAWS Cruiser - 16inch - Auto Mounting - Standard Hand Control
WCA872130	PAWS Cruiser - 16inch - Auto Mounting - Quad/Tetra Hand Control